[Gebruiksaanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Transferpette® S

Pipettierhelfer | Pipette controllers

Impressum

BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

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Table of contents

1 Using the operating manual	35		
1.1 Signal words and their meaning	35		
1.2 Presentation of descriptions of actions	35		
2 Safety regulations.....	36		
2.1 Safety regulations	36		
2.2 Target group.....	37		
2.3 Use	37		
2.4 Purpose	37		
2.5 Improper use	37		
2.6 Foreseeable misuse	37		
2.7 Limitations of Use	37		
2.8 Application restrictions	37		
2.9 Operating exclusions	38		
3 Scope of delivery	39		
4 Functional and operational components	40		
5 Pipetting	41		
6 Checking the volume	44		
7 Accuracy table	45		
8 Adjustment – Easy Calibration.....	47		
9 Disinfection/Autoclaving	49		
9.1 Autoclaving.....	49		
9.2 UV disinfection	49		
9.3 PE filter	49		
10 Maintenance.....	50		
10.1 Check pipetting unit	50		
10.2 Disassembling/cleaning (up to 1,000 µl)	51		
10.3 Disassembly/cleaning (2–10 ml)	52		
11 Troubleshooting	53		
12 Product markings	54		
13 Order Information.....	55		
14 Repairs.....	59		
		14.1 Sending for repair	59
15 Calibration service	61		
16 Information about your laboratory device.....	62		
17 Warranty.....	63		
18 Disposal	64		

1 Using the operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the instructions for use at www.brand.de/om.

1.1 Signal words and their meaning

 WARNING	The signal word WARNING indicates a hazardous situation that, if not avoided, could result in death or serious injury.
or  WARNING! ...	
 CAUTION	The signal word CAUTION indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
or  CAUTION! ...	
NOTICE	The signal word NOTICE is used to address practices not related to physical injury. Example: Possible property damage.
or NOTICE ...	shows operation tips or other helpful information.

1.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 Safety regulations

Please read carefully!

The Transferpette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the operating manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations (e.g. wear protective clothing, eye protection and protective gloves).
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/media must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.
- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the "Troubleshooting" section in the operating manual. Contact the manufacturer, if necessary.
- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the operating manual!
- Only autoclave the laboratory device when this is permissible according to the operating manual.

2.2 Target group

The operating manual is intended for users who use the laboratory instrument in the course of their professional activities. Users are familiar with the typical safety regulations and working methods in laboratories and have been trained accordingly. They can recognize potential hazards and protect themselves from them. The operating manual assumes this expertise and does not replace basic laboratory training or specific safety training.

2.3 Use

Use the laboratory device Transferpette® S only for pipetting or dispensing liquids within the defined limits of use.

2.4 Purpose

This is an air displacement pipette for pipetting liquids of low to medium viscosity.

2.5 Improper use

Various risks may arise if the laboratory instrument is used improperly. These risks include: inaccurate liquid delivery, damage to the laboratory instrument, and the risk of contamination, infection, and injury from contact with the pipetted media.

Any use other than for pipetting or dispensing liquids within the defined operating limits is considered improper use.

2.6 Foreseeable misuse

A typical misuse is pipetting or dispensing liquids with too high viscosity or using unsuitable tips. Furthermore, setting the volume while the volume-change protection is locked is a foreseeable misuse.

2.7 Limitations of Use

The pipette is used for dispensing liquids within the following limits:

- Operating temperature from +15 °C to +40 °C (59 °F to 104 °F). Additional temperatures upon request
- vapor pressure up to 500 mbar
- Viscosity: 260 mPa s

With viscous liquids, the dispensing speed may need to be adjusted.

2.8 Application restrictions

Viscous and wetting liquids may compromise volumetric accuracy. Volumetric accuracy may also be affected when pipetting liquids whose temperature deviates from the ambient temperature by more than ± 1 °C/ ± 1.8 °F.

2.9 Operating exclusions

The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:


- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Fluoroelastomer rubber (FKM)
 - Polycarbonate (inspection window)
 - Polyetheretherketone (PEEK)
 - Polyoxymethylene (POM)
 - Polyphenylene sulfide (PPS) (with variable 50 µl instrument, 5 ml, 10 ml pipetting unit)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)

Additional information on the chemical resistance of plastics can be found www.brand.de

3 Scope of delivery

Transferpette® S, with quality certificate and silicone grease.

4 Functional and operational components

- | | |
|------------------------------------|-------------------------------|
| 1 Pipetting button | 2 Volume-setting wheel |
| 3 Easy Calibration function | 4 Finger rest |
| 5 Pipetting shaft | 6 Tip cone |
| 7 Handle | 8 Volume display |
| 9 Volume-change protection | 10 Tip ejection key |

Label window

The instrument can be individually labeled on the finger rest:

- a. Remove the label window on the finger rest.
- b. Mark the labeling film.
- c. Reinsert the labeling film with window.


5 Pipetting

1. Attach tip

2 ml, 5 ml and 10 ml devices should only be used with a built-in PE filter; see UV disinfection, p. 49.


Perfect analysis results can only be achieved by using quality tips. We recommend BRAND pipette tips. For additional information, refer to the accuracy table Accuracy table, p. 45.

Pipette tips are disposable products!


- a. Insert tips vertically:
Use the correct tips, in accordance with the volume range or color code!
Make sure that the tips are firmly in place and leak tight.

2. Set volume


- a. Slide the volume-change protection upward (UNLOCK).
- b. Turn the volume-setting wheel to select the desired volume. In doing so, turn the adjustment wheel steadily, avoiding abrupt turning motions.


- c. Slide the volume-change protection downward (LOCK). The volume-setting wheel becomes noticeably more difficult to turn, but movement is not completely blocked. Avoid setting the volume while the volume-change protection is locked. This is a form of misuse that can lead to volume deviations. If this occurs, perform a volume check and, if necessary, adjust the instrument using EasyCalibration.

3. Take sample

We recommend rinsing the pipette tip 5 times with the liquid before pipetting (aspirate and dispense liquid 5 times) to achieve maximum precision and accuracy. See also Checking the volume, p. 44.

- Press the pipetting button until first resistance is felt.
- Hold the instrument vertically and immerse the tip in the liquid.


- Allow the pipetting button to steadily move back to its original position.

Leave the tip immersed in the liquid for a few seconds, so that the set volume is aspirated completely. This is especially important when pipetting viscous media and when using pipettes with large volumes.


NOTICE Do not lay the device down when the tip is filled; this can cause the medium to flow into the device and contaminate it!

Volume range	Immersion depth [mm]	Wait time [s]
0.1–1 µl	1–2	1
1–100 µl	2–3	1
100–1000 µl	2–4	1
> 1,000 µl	3–6	3

4. Dispense sample


- a. Place the pipette tip against the vessel wall. Hold the pipette at an angle of 30–45° to the vessel wall.
- b. Press the pipetting button at a uniform speed until the first resistance is felt and hold it. To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density media.

- c. Empty the tip completely by overtravel: Press the pipetting button to the second stop.
- d. While doing this, wipe the pipette tip against the vessel wall.
- e. Remove the pipette tip from the vessel wall and allow the pipetting button to move back to its original position.

5. Discard tip

NOTICE The instrument should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

- a. Hold the pipette shaft over a suitable disposal bin and press the tip ejection button all the way down.

6 Checking the volume

We recommend testing the device every 3 to 12 months, depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de.

Gravimetric volume testing of the pipette is carried out according to the following steps and complies with ISO 8655.

- Set the maximum specified instrument volume (for procedure, see Pipetting, p. 41).
- Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.
- Aspirate the test liquid and pipette into the weighing vessel.
- Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- At least 10 pipetting series and weighings in three volume ranges (100%, 50%, 10%) are recommended. Two tips must be used for each volume range to be tested.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

7 Accuracy table

Transferpette® S, adjustable

Volume range [μl]	Partial volume [μl]	A* ≤ ± %	CV* ≤ %	Sub steps [μl]	Recommended tip type [μl]
0.1 - 1	1 0.5 0.1	2 4 20	1.2 2.4 12	0.001	0.1 - 20
0.1 - 2.5	2.5 1.25 0.25	1.4 2.5 12	0.7 1.5 6	0.002	0.5 - 20
0.5 - 10	10 5 1	1 1.6 7	0.5 1 4	0.01	0.5 - 20
2 - 20	20 10 2	0.8 1.2 5	0.4 0.7 2	0.02	2 - 200
5 - 50	50 25 5	0.8 1.2 4	0.3 0.5 2	0.05	2 - 200
10 - 100	100 50 10	0.6 0.8 3	0.2 0.4 1	0.1	2 - 200
20 - 200	200 100 20	0.6 0.8 3	0.2 0.3 0.6	0.2	2 - 200
100 - 1000	1000 500 100	0.6 0.8 3	0.2 0.3 0.6	1	50 - 1000
500...5,000	5000 2500 500	0.6 0.8 3	0.2 0.3 0.6	5	500...5,000
1000 - 10,000	10000 5000 1000	0.6 0.8 3	0.2 0.3 0.6	10	1,000...10000

*A = Accuracy, CV = Coefficient of Variation

Transferpette® S, fixed-volume

Volume range [μl]	A* ≤ ± %	CV* ≤ %	Recommended tip type [μl]
10	1	0.5	0.5 - 20
20	0.8	0.4	2 - 200
25	0.8	0.4	2 - 200
50	0.8	0.4	2 - 200
100	0.6	0.2	2 - 200
200	0.6	0.2	2 - 200
500	0.6	0.2	50 - 1000
1000	0.6	0.2	50 - 1000
2000	0.8	0.3	500...5,000


*A = Accuracy, CV = Coefficient of Variation


Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

8 Adjustment – Easy Calibration

The device is permanently calibrated for aqueous solutions. If it is determined that the pipette is operating inaccurately or to adjust the device to work with solutions of varying density and viscosity or with specially-shaped pipette tips, it can be calibrated using the Easy Calibration Technique.

- a. Perform a volume check and determine the actual value; see Checking the volume.
- b. Remove the label window (1) and labeling film (2): lift the hook slightly and remove it upward.
- c. Using a paper clip or an unused pipette tip, remove the protective film (3) (the protective film can be discarded).

- d. Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

- e. Set the adjustment volume:
Transferpette® S, adjustable: with the volume-setting wheel in the UNLOCK position, set to the previously determined actual value.
Transferpette® S, fixed-volume: set the volume by rotating in the +/- direction.
 A volume check is recommended after every adjustment.

f.

- f. Slide the adjustment slider completely back again, push the volume-setting wheel downward and release the adjustment slider. Re-attach the labeling film and re-assemble the label window.

NOTICE The change to factory settings is indicated by the red adjustment slider now visible in the label window.

9 Disinfection/Autoclaving

9.1 Autoclaving

NOTICE Check the effectiveness of autoclaving yourself!

Maximum safety is achieved through vacuum sterilization. We recommend using sterilization bags.

- Eject the pipette tip.
- Fully unlock the volume adjustment lock (UNLOCK).
- Set the volume to the full value (e.g. 11.25 or 11.26) and leave the volume adjustment lock unlocked. This prevents jamming or damaging the volume setting during autoclaving.

Correct

Incorrect

- Pack the instrument in a sterilization bag, observing any applicable packaging regulations.
- Autoclave the complete pipette without further disassembly. Recommendation for autoclaving according to DIN EN 285 (see table below).
- Allow the pipette to completely cool and dry.

Temperature:	121 °C
Pressure	2 bar
Holding time in autoclaves	15 min

After autoclaving, tighten the connection between the hand grip and the pipette shaft if necessary.

9.2 UV disinfection

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

9.3 PE filter

PE filter for Transferpette® S 2 ml, 5 ml + 10 ml:

A hydrophobic PE filter is used to prevent liquid from entering the pipette.

Change the filter if it becomes wet or contaminated.

- Use a flat object, such as a screwdriver.
- Remove the filter carefully, without damaging the tip cone.

Remove the filter before autoclaving!

The device can also be operated without a filter.

10 Maintenance

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.


We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see --- FEHLENDER LINK ---.

10.1 Check pipetting unit

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.

We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see Troubleshooting, p. 53.

10.2 Disassembling/cleaning (up to 1,000 µl)


- a. Detach the pipette shaft (S) from the hand grip by unscrewing it.
- b. Unscrew the upper part of the ejector unit (A) from the pipette shaft.
- c. Pull out the shaft (B, C and D) from the lower part of the ejector unit (E).
- d. Unscrew the piston unit (B).

NOTICE Do not disassemble the piston unit (B) any further!

- e. Remove the seal with spring (C) (not possible on Transferpette® S 1 µl, 2.5 µl, and 10 µl!).
- f. Clean the parts shown with a soap solution or isopropanol, and then rinse with distilled water.
- g. Dry the parts (max. 120 °C/248 °F).
- h. Grease piston and seal with a very thin layer of supplied silicone grease.

Reassemble the cooled parts in reverse order. Only hand-tighten the piston unit and the upper part of the ejector unit (A, B).

10.3 Disassembly/cleaning (2–10 ml)

- Remove the entire shaft (S) from the hand grip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
- Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
- Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
- Remove the O-ring-seal from the piston unit and clean it.

NOTICE Do not disassemble the piston unit (G) any further!

- Clean the piston unit (G) and the bottom part of the shaft (H) with a soap solution or isopropanol, and then rinse with distilled water.
- Dry the parts (max. 120 °C/248 °F) and allow them to cool.
- Carefully lubricate the inside and outside of the O-ring (O) and mount it on the piston.

Reassemble the individual components in reverse order.

11 Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (device leaking)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Firmly press tip on
The instrument does not aspirate or aspirates too little; the dispensed volume is too low	Seal contaminated	Clean seal
	Seal or cone is damaged	Replace seal or shaft
	Piston is contaminated or damaged	Clean or replace piston
Aspiration is very slow	Shaft is clogged	Clean shaft
	Filter contaminated on 2 ml, 5 ml or 10 ml devices	Change filter
Dispensed volume too large	Pipetting button pressed too far (to the over-stroke point) before aspirating	Ensure proper handling.
Piston sluggish	Piston is contaminated or not greased	Clean piston and apply grease

12 Product markings

Symbol or number	Meaning
	Read the user manual.
XXZXXXXX	Serial number
 25	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
Data matrix code or quick response code	The codes link to the BRAND MyProduct website.
www.brand.de/ip	Hyperlink to the BRAND patent page

13 Order Information

Transferpette® S, fixed-volume


Volume	Item no.
10 µl	705808
20 µl	705816
25 µl	705820
50 µl	705828
100 µl	705838
200 µl	705844
500 µl	705854
1,000 µl	705862
2000 µl	705864

Transferpette® S, adjustable

Volume	Item no.
0.1–1 µl	705868
0.1–2.5 µl	705869
0.5–10 µl	705870
2–20 µl	705872
5–50 µl	705873
10–100 µl	705874
20–200 µl	705878
100–1,000 µl	705880
0.5–5 ml	705882
1–10 ml	705884

Spare parts - Transferpette® S up to 1,000 µl

Transferpette® S, fixed-volume


Volume	A	B	C	D	E
10 µl	705508	704601	–	704721 *	704739
20 µl	705509	704602	704610	704723	704740
25 µl	705509	704608	704614	704723	704741
50 µl	705509	704654	704661	704724	704742
100 µl	705509	704654	704661	704724	704743
200 µl	705509	704655	704662	704725	704745
500 µl	705511	704656	704663	704726	704746
1,000 µl	705511	704656	704663	704726	704747

Transferpette® S, adjustable

Volume	A	B	C	D	E
0.1– 1 µl	705508	704600	–	704718 *	704730
0.1– 2.5 µl	705508	704616	–	704719 *	704731
0.5– 10 µl	705508	704601	–	704721 *	704732
2– 20 µl	705509	704602	704610	704723	704733
5– 50 µl	705509	704615	704617	704722	704734
10– 100 µl	705509	704654	704661	704724	704735
20– 200 µl	705509	704655	704662	704725	704736
100– 1000 µl	705511	704656	704663	704726	704737

* Seal permanently installed in shaft – not removable!

Spare parts - Transferpette® S, 2 ml, 5 ml, and 10 ml

NOTICE The appearance and dimensions of the spare parts correspond to the respective nominal volume.

Transferpette® S, fixed-volume

Volume	F + F'	G	H	I	O
2 ml	704765	704606	703247	704626	7288
2 ml (from 01/2026)	704765	704669	704689	704626	7228

Transferpette® S, adjustable

Volume	F + F'	G	H	I	O
0.5–5 ml	704766	704606	703247	704626	7288
1–10 ml	704767	704607	704628	704626	7298

Accessories

Table stand for 1 pipette
Item no. [703440](#)

Wall mount
Item no. [704882](#)

Table stand for 6 single-channel or multi-channel pipettes
Item no. [704807](#)

Shelf/rack mount
Item no. [704881](#)

Label window
Item no. [704750](#)

Labeling film
Item no. [704751](#)

Filter for volume range 2–5 ml, pkg. unit 25 pc
Item no. [704652](#)

Silicone grease
volume range up to 1000 µl
Item no. [705502](#)

Silicone grease
for volumes 2 ml/5 ml/
10 ml
Item no. [703677](#)

Filter for volume range up to 10 ml, pkg. unit 25 pc
Item no. [704653](#)

PLT unit (pipette leak detector)
Item no. [703970](#)

14 Repairs

14.1 Sending for repair

NOTICE Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1-860-767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

15 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the instrument. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive liquids.

The detailed testing instructions are available for download on www.brand.de or www.brandtech.com.

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or our accredited calibration laboratory. Simply send us the instruments to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days. The instruments are supplied with a detailed calibration certificate or a calibration slip in accordance with DIN EN ISO/IEC 17025. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on www.brand.de (see Service and Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. If a factory calibration is requested, the service partner can forward the instrument to BRAND.

16 Information about your laboratory device

The online service MyProduct (<https://www.brand.de/myproduct>) offers quality certificates, accessories, and technical documentation for your laboratory instrument Transferpette® S. By entering the serial and item number there, you will receive information specific to your individual instrument.

On Transferpette® S, you will also find serialized information encoded in the quick response code. This also links to the MyProduct website, where you can access this operating manual and additional certificates for your instrument.

You will also find a data matrix code on some instruments (Transferpette® S, HandyStep® touch, and HandyStep touch® S). Scan this using a standard reader app to access the above-mentioned information via the URL <https://www.brand.de/myproduct>.

17 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

18 Disposal

At the end of their service life, ensure the instruments are properly disposed of and, where possible, recycled. Disposal must be carried out in accordance with applicable national and local regulations.

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